

Oral Rehabilitation Device

Littoremeter



These items are for those with deteriorated oral function due to aging, disease or accident. 3 types of devices are available. These items help to improve the oral function by training oral muscle and tongue.



* Orbicular muscle of the mouth can be trained with this device by placing the mouth piece between the front teeth and lips and pulling it. The measure on the device provides a numerical guideline of muscle strength.



* The labial closing force can be trained with this device by placing upper and lower lips onto the white guide and repeatedly tightening it.
* Three levels of load can be chosen by changing the inner spring (100g/200g/400g)

Agero M Medical



* These devices are for training lingual muscle by placing the spoon shape mouth piece on the tongue and repeatedly lifting it up.
* Two sizes are available (M and S size)
* Three levels of load can be chosen by changing inner spring (100g/200g/400g)

